



## PharmaScoops™

kersting PharmaScoops™ werden aus hochwertigem V4A Edelstahl hergestellt.

Durch die gratfreie Konstruktion sind PharmaScoops™ für Anwendungen mit höchsten Hygiene-Anforderungen geeignet.



Linke Seite: PharmaScoops™ Übersicht

Recht Seite: PharmaScoop™ rechteckig (278.1500)

**PharmaScoops™ sind sehr robust und können bei guter Pflege viele Jahre im Einsatz bleiben.**

Die Reinigungsmethode hängt von der Art der Produkte ab, mit denen die PharmaScoops™ in Berührung kommen. Bitte stellen Sie sicher, dass sie nach der Reinigung gründlich getrocknet werden.

**PharmaScoops™ können auf folgende Arten gereinigt werden:**

- Spülen per Hand
- Spülmaschine
- Ultraschallbad
- Autoklave 121°C

**Bitte reinigen Sie die PharmaScoop™ vor dem ersten Gebrauch.**

Artikel Nr.	Beschreibung	Volumen	Ø mm	Gesamtlänge mm	Schaufellänge mm
279.0010	Pharma Scoop	10 ml	25	75	45
279.0050	Pharma Scoop	50 ml	38	110	70
279.0100	Pharma Scoop	100 ml	50	180	100
279.0500	Pharma Scoop	500 ml	95	265	140
279.0750	Pharma Scoop	750 ml	95	330	210
279.1100	Pharma Scoop	1000 ml	110	340	235
279.1200	Pharma Scoop	1200 ml	125	305	195
279.1700	Pharma Scoop	1700 ml	135	355	230
279.2500	Pharma Scoop	2500 ml	185	320	215
278.1500	Pharma Scoop rechteckig	1500 ml	160	310	190

### KERSTING GMBH SAMPLING + GROUNDING

Gallbergweg 21, 59929 Brilon / Germany


☎ + 49 2961 7405 – 302, F: + 49 2961 7405 – 349, mailto: [info@kersting-ind.de](mailto:info@kersting-ind.de), [www.kersting-ind.de](http://www.kersting-ind.de)



## **PharmaScoops™**

kersting *PharmaScoops™* are made from high quality 316 stainless steel.

Their crevice free construction makes them perfectly suitable for applications in hygiene conscious industries.



Left side: PharmaScoops™ range overview

Right side: PharmaScoop™ rectangular (278.1500)

*PharmaScoops™* are very robust and with proper care will last many years.

The cleaning method of choice will depend upon the nature of the products with which the *PharmaScoop™* has been used. Ensure that the *PharmaScoop™* is thoroughly dried after washing.

The following cleaning methods can be used for the *PharmaScoops™*:

- Hand Wash
- Washing Machine
- Ultrasonic Bath
- Autoclave 121°C

Please clean the *PharmaScoop™* before it is used for the first time.

Item no.	Description	Volume	Ø mm	Total Length mm	Scoop Length mm
279.0010	Pharma Scoop	10 ml	25	75	45
279.0050	Pharma Scoop	50 ml	38	110	70
279.0100	Pharma Scoop	100 ml	50	180	100
279.0500	Pharma Scoop	500 ml	95	265	140
279.0750	Pharma Scoop	750 ml	95	330	210
279.1100	Pharma Scoop	1000 ml	110	340	235
279.1200	Pharma Scoop	1200 ml	125	305	195
279.1700	Pharma Scoop	1700 ml	135	355	230
279.2500	Pharma Scoop	2500 ml	185	320	215
278.1500	Pharma Scoop rectangular	1500 ml	160	310	190

**KERSTING GMBH SAMPLING + GROUNDING**

Gallbergweg 21, 59929 Brilon / Germany

☎ + 49 2961 7405 – 302, F: + 49 2961 7405 – 349, mailto: [info@kersting-ind.de](mailto:info@kersting-ind.de), [www.kersting-ind.de](http://www.kersting-ind.de)